



AMERICAN INTERNATIONAL UNIVERSITY – BANGLADESH

Department of Natural Science (Chemistry)
Faculty of Science & Technology
Programs: B.Sc. Eng'g (EEE/CSE/IPE)
CHEM 1101: CHEMISTRY

Chemistry Lab Report

Semester: Spring

Session: 2022-2023

NO EXPERIMENT, NO REPORT

Experiment No: 2 Standardization of Hydrochloric Acid (HCl)
Name of the Experiment: Standardization of Hydrochloric Acid (HCl)
Solution with Standard Sodium Hydroxide (NaOH) Solution.

Date of Performance: 15-02-23, Date of Submission: 22-02-23

Course-Teacher: Dr. Mohammad Tariqul Islam

Instructions:

1. A lab report consists of three parts: a cover page, body of the report and a data and results sheet (lab-sheet).
2. This is the cover page of a report and students will collect and preserve the lab-sheet of a particular experiment to be performed.
3. Body of the report includes-(1) Objective of the Experiment, (2) Theory, (3) Name of the Chemicals, (4) Name of the Apparatus, (5) Percentage of Error (if necessary) and (6) Discussion (I. Precautions taken, II. Possible errors).
4. Use A₄-size off-set paper, write on one side of the paper by hand keeping suitable margin.
5. Staple the lab-sheet at the end of the report and cover page on the top.
6. Submit the report in time to avoid deduction of marks.
7. Students working in a group will write and submit the report individually.
8. Copying of the report from others is strictly prohibited.

Name of the Student: TRIDIB SARKAR
ID No: 22-46449-1, Section: F, Group: 3

FOR FACULTY USE ONLY


Faculty comments: Signature:
Date:

Objective: To know the strength of HCl solution (being a solution made from secondary standard substance) against a previously standard solution by acid-base titration.


Theory:

(i) Method: In presence of a suitable indicator the volumetric analysis in which a standard solution is added in another solution (whose strength is not known) to reach its end point to determine the strength of that solution is called titration. A solution of known concentration is called a standard solution. In titration, primary and secondary standard solutions are commonly used.

(ii) Reactions:




This is an acid-base neutralization reaction. Oxalic acid reacted with NaOH and produce $\text{NaO}_2\text{C}-\text{CO}_2\text{Na}$ and water.




This is also an acid-base neutralization reaction. Here sodium hydroxide reacts with hydrochloric acid and produce sodium chloride (NaCl) and water (H_2O).

(iii) Indicators:

In first reaction $\text{H}_2\text{C}-\text{CO}_2\text{H}$ is weak acid and NaOH is strong base. So, here used phenolphthalein. In acid it is colourless and in base it is pink.



In second reaction NaOH is strong base and HCl is strong acid so any indicator is preferable. Here in this reaction methyl orange indicator is used.



Required chemicals:

<u>Name of the chemicals</u>	<u>Chemical formula</u>
1. Supplied NaOH solution	NaOH
2. Standard Oxalic Acid Solution	C ₂ H ₂ O ₄
3. Phenolphthalein indicator	C ₂₀ H ₁₄ O ₄
4. HCl acid solution	HCl
5. Methyl orange indicator	C ₁₄ H ₁₄ N ₃ NaO ₃ S

Apparatus:

- | | |
|--------------------------------|-----------------------------|
| (i) Burette (50 ml) | (v) Watch glass. |
| (ii) Pipette (10 ml) | (vi) Pipette filler. |
| (iii) Conical flask (250 ml) | (vii) Dropper. |
| (iv) Volumetric flask (100 ml) | (viii) Stand and clamp etc. |

Name: TRIDIB SARKAR ID No: 22-A/AAA-1, Section (Group): F(3)

AMERICAN
INTERNATIONAL
UNIVERSITY –
BANGLADESH (AIUB)



Experiment 2

CHEM 1101: CHEMISTRY (EEE/CoE/CSE/IPE)

EXPERIMENT NO. 2: STANDARDIZATION OF HYDROCHLORIC ACID (HCl) SOLUTION WITH STANDARD SODIUM HYDROXIDE (NaOH) SOLUTION.

OBJECTIVE: To know the strength of HCl solution (being a solution made from secondary standard substance) against a previously standard solution by acid-base titration.

THEORY:

- (i) *Method:* Acid-base titration
- (ii) *Reactions:* 1. $\text{HO}_2\text{C}-\text{CO}_2\text{H} + 2\text{NaOH} = \text{NaO}_2\text{C}-\text{CO}_2\text{Na} + 2\text{H}_2\text{O}$
2. $\text{NaOH} + \text{HCl} = \text{NaCl} + \text{H}_2\text{O}$
- (iii) *Indicators:* Phenolphthalein, Methyl orange

APPARATUS:

Burette (50mL), pipette (10mL), conical flask (250mL), volumetric flask (100mL), watch glass, pipette filler, dropper, Stand and clamp etc.

REQUIRED CHEMICALS:

1. Supplied NaOH solution
2. Standard oxalic acid solution
3. HCl acid solution
4. Phenolphthalein indicator
5. Methyl orange indicator

(A) Standardize the supplied NaOH solution as in Experiment No. 1

$$\text{Strength of oxalic acid solution} = \frac{\text{Weight taken (in gm)} \times 0.1}{0.63} (\text{N}) = 0.10158 \text{ N}$$

Table-1: Standardization of supplied NaOH solution against standard oxalic acid solution by acid-base titration.

No. of reading	Vol. of NaOH (in mL)	Vol. of Oxalic acid (burette reading) (in mL)			Mean (in mL)
		Initial	Final	Difference	
1	10	0	9.40	9.40	
2	10	9.4	18.60	9.20	
3	10	18.6			9.30

Strength of supplied NaOH solution:

$$V_{\text{NaOH}} \times N_{\text{NaOH}} = V_{\text{Oxalic acid}} \times N_{\text{Oxalic acid}}$$

$$\Rightarrow 10 \times N_{\text{NaOH}} = 9.30 \times 0.10158$$

$$\therefore N_{\text{NaOH}} = 0.0943$$

Chemistry Lab Sheet

Name: TRIDIB SARKAR..... ID No. 22-16444-1, Section (Group): F (3)

(Expt. 2 contd.)

(B) Preparation of approximately 0.1N hydrochloric acid solution:

Take 10 ml conc. HCl in a 1000 ml measuring flask and add distilled water up to the mark.

PROCEDURE: Take 10 mL of NaOH solution in a conical flask by means of a pipette and dilute it to about 50 mL. Add 2-3 drops of methyl orange indicator to the solution. Then add previously prepared (approx. 0.1N) HCl acid solution drop wise from a burette. Shake the flask frequently during addition of HCl acid. Stop the addition of HCl acid solution as soon as the yellow color of the solution just changes to orange or pink. Note the burette reading. Repeat the process at least three times and take the mean of the readings. Calculate the strength of the dilute HCl solution and from there calculate the strength of commercial HCl.

EXPERIMENTAL DATA:

Table-2: Standardization of supplied HCl solution against standard NaOH solution by acid-base titration.

No. of reading	Vol. of NaOH (in mL)	Vol. of HCl (burette reading) (in mL)			Mean (in mL)
		Initial	Final	Difference	
1	10	17.40	28.40	11.00	
2	10	28.40	38.80	10.40	
3	10	20.90	32.50	11.60	
4	10				11.00

CALCULATIONS:

(A) Strength of supplied dil. HCl solution:

$$V_{\text{NaOH}} \times N_{\text{NaOH}} = V_{\text{dil. HCl}} \times N_{\text{dil. HCl}} \text{ to be determined}$$

$$10 \times 0.0943 = 11.00 \times N$$

$$\Rightarrow N = 0.08575 \Rightarrow 0.09$$

(B) Strength of conc. HCl solution:

$$V_{\text{dil. HCl}} \times N_{\text{dil. HCl determined}} = V_{\text{conc. HCl taken}} \times N_{\text{conc. HCl to be determined}}$$

$$1000 \times 0.08575 = 10 \times N$$

$$\Rightarrow N = 8.5750$$

RESULTS:

(a) The strength of supplied dil. HCl solution is 0.09 N

(b) The strength of conc. HCl solution is 8.575 N

- What is normality and molarity?
- Atomic weight, molecular weight and gram equivalent weight of NaOH, HCl and HOOC-COOH, $2\text{H}_2\text{O}$
- Why phenolphthalein and/or methyl orange are used?
- Reason of using methyl orange instead of phenolphthalein.

Text: M. Mahbubul Huque and A. Jabber Mian, "Practical Chemistry", 2nd ed. (1972)

Discussion:

Precautions:

1. All the apparatus were washed with first normal water and then with distilled water.
2. Had to be very keen while adding oxalic acid solution and HCl acid solution.
3. Solution was shaken spontaneously.

Possible errors:

1. Error might be occurred while taking the burette reading.
2. Might added more than 40 ml distilke water in the NaOH solution.
3. The lower meniscus and the eyelevel should be at same level while taken the burette reading.